

IHE Integration Statement			
Vendor	Product Name	Version	Date
lifeIMAGE	lifeIMAGE	4.0	1/29/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Consistent Time	Time Client	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V2)
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V3)
Cross-Enterprise Document Sharing	Document Repository	None
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing for Imaging	Document Registry	None
Cross-Enterprise Document Sharing for Imaging	Document Repository	None

Internet address for vendors IHE information :<http://www.lifeimage.com>

Links to Standards Conformance Statements for the Implementation

**HL7**

**DICOM**

Link to general information on IHE :[www.ihe.net](http://www.ihe.net)